


STAGE MICROMETER - S-14770



Item no: S-14770 - Stage Micrometer, 1mm total scale with 200 divisions, 0.005mm/division

Glass: 1.0" x 3.0" x .110" thick
Lines: Chrome deposition
Flatness: 2 Waves or better
Tolerances: $\pm 0.005\text{mm}$ unless otherwise noted
Centering: Within 0.005"
Surface Quality: 20/10
Parallelism: 2 min.

Calibration traceable to NIST available at additional cost

Colusa Science

PH: 800-830-1885 • FAX: 800-830-1887 • EMAIL: info@col-sci.com
327 College St., Ste 203 • Woodland, CA. 95695 USA