





Ph: 888-777-8835 • 209-839-2990 • Fax: 209-839-2999  
sales@m-depot.com • www.microscope-depot.com  
392 W. Larch Rd., Bldg. 28, Tracy, CA. 95304 USA



KR-806 / 1" Stage Micrometer / .1" - .01" - .001" Div. / Glass Slide -- 1" x 3" x .125" thick



KR-812 / 25mm Stage Micrometer / 1.0mm - .1mm - .02mm Div. / Glass Slide -- 1" x 3" x .125" thick

KR-814 / 1" & 25mm Stage Micrometer / both above scales on same slide / Glass Slide -- 1" x 3" x .125" thick